Color Video Printer

Instructions For Use Page 2

Before operating this unit, please read this manual thoroughly and retain it for future reference.

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UP-1850EPM

WARNING

To prevent fire or shock hazard, do not expose the unit to rain or moisture.

To avoid electrical shock, do not open the cabinet. Refer servicing to qualified personnel only.

Symbols on the products

This symbol indicates type B equipment classified in accordance with IEC Publication 601-1 Safety of medical electrical equipment.


This symbol indicates the equipotential terminal which brings the various parts of a system to the same potential.

Important safeguards/notices for use in the medical environments

1. All the equipments connected to this unit shall be certified according to Standard IEC601-1, IEC950, IEC65 or other IEC/ISO Standards applicable to the equipments.
2. When this unit is used together with other equipment in the patient area*, the equipment shall be either powered by an isolation transformer or connected via an additional protective earth terminal to system ground unless it is certified according to Standard IEC601-1.

* Patient Area

3. The leakage current could increase when connected to other equipment.

For the customers in the U.S.A.

Warning

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

You are cautioned that any changes or modifications not expressly approved in this manual could void your authority to operate this equipment.

The shielded interface cable recommended in this manual must be used with this equipment in order to comply with the limits for a digital device pursuant to Subpart B of Part 15 of FCC Rules.

For the customers in the United Kingdom

WARNING

THIS APPARATUS MUST BE EARTHED

IMPORTANT

The wires in this mains lead are coloured in accordance with the following code:

Green-and-yellow:	Earth
Blue:	Neutral
Brown:	Live

As the colours of the wires in the mains lead of this apparatus may not correspond with the coloured markings identifying the terminals in your plug proceed as follows:

The wire which is coloured green-and-yellow must be connected to the terminal in the plug which is marked by the letter E or by the safety earth symbol \perp or coloured green or green-and-yellow.

The wire which is coloured blue must be connected to the terminal which is marked with the letter N or coloured black.

The wire which is coloured brown must be connected to the terminal which is marked with the letter L or coloured red.

This PRINTER is designed to use CARTRIDGES bearing the VP trademark to yield quality prints.

VP and Design is a trademark used under license.

Table of Contents

Introduction

About This Manual 5

System Overview 6

System Configuration 7

Operation

Before Printing 8

Loading an Ink Ribbon Cassette 8

Loading Paper 10

Selecting the Input Signal 12

Making Full-Size Printouts 14

Making Multiple Copies of Identical Image 18

Making Variations of Printouts 20

Printing Operation Flowchart 21

About Memory 22

Selecting the Memory Mode 26

Selecting the Print Mode 28

Making Printouts of Multiple Reduced Images 31

Making Printouts of Multiple Images in Multiple Memory Pages 35

Making Printouts with a Caption 40

Entering a Caption 41

Deleting the Images Stored into Memory Pages 45

Simultaneously Deleting All Images in a Single Memory Page 45

Deleting Part of Multiple Images in a Single Memory Page 48

Installation and Adjustment

Supplied Accessories 51

Connections 52

Making Connections for Storing Video Images 52

Making Connections for Viewing Images to be Printed on the Video Monitor 53

Making Connections to Enable Remote Control 54

Preparing the Remote Commanders 55

Using the Supplied Remote Commander 55

Using the Remote Commander (Not Supplied) 56

Adjusting the Printout Quality 57

Adjusting the Sharpness 57

Adjusting the Printout Color 58

Table of Contents (continued)

Installation and Adjustment

Printer Initial Setup 63

- Changing the Printout Area 63
- Compensating for the Input Signals 66
- Matching the Video Monitor Color to the Printer Color 69
- Erasing the Screen Display 75
- Selecting Operation Mode for Automatic Printing Capabilities 78
- Setting the Dip Switches 84

Others

Precautions 88

- Safety 88
- Installation 88
- Cleaning 89

Ink Ribbon Cassette and Paper 90

Specifications 91

Error Messages 93

- If the Paper Jams 94

Troubleshooting 97

Location and Function of Parts and

Controls 98

- Front 98
- Rear 100
- Remote Commander RM-5100 101
- Comparing Printer Window Display with the Video Monitor Display 102
- Monitor Display 103

Index 107

About This Manual

This manual covers the following UP-1800EPM series color video printers.

UP-1800EPM (single-frame memory)

UP-1850EPM (three-frame memories)

Wherever the operation or any other item differs between the models, this manual clearly describes those differences.

Organization of this manual

Introduction

Describes the features and system configuration of the color video printer.

Operation

Describes actual printing operation once all connections and adjustments have been made, as explained in the next chapter. You will be able to make various types of printouts after reading through this chapter.

Installation and adjustment

Describes how to make connections and make adjustments using the menus displayed on the video monitor and printer window display. Once all connections and adjustments have been made, there should be no need to perform these operations again during normal printing operations. These operations must, however, be performed after reinstalling, or if the picture quality degrades, or if adjustment becomes necessary because the peripheral equipments is changed. Also, covered is the use of the printer's remote commander.

Others

Notes the precautions to be observed when using the printer, lists errors and their handling, and explains troubleshooting. Also provided is information on the locations and functions of parts and controls, and the on-screen messages and menus used to operate the printer. Should you encounter any unfamiliar terms or items while reading this manual, consult the index at the end of the manual.

Conventions used

Cross references

Throughout this manual you will find the references to other sections of the manual that contain related information.

Important note

Be sure to read the sections of the manual marked **Note**. They explain points that you should be aware of to operate the printer correctly and prevent malfunctions.

Index


Use the index, in addition to the table of contents, to find information you need when using the printer.

System Overview


The Sony UP-1800EPM/1850EPM color video printer is designed for capturing images from video equipment and making printouts of those images. By changing the printout mode, different types of printouts can be made. Also, you can add a caption to a printout. Printer setup is done interactively by picking from displayed menus and is also possible by DIP switch setting. The printer can make the following types of printouts.

Printouts that can be made with the printer


Printout of a full-size image (page 14)


Printouts of mirror images (page 28)


Printout of four reduced images (pages 26 and 31)


Printouts with an insert (page 33)

Printout of 16 reduced images (pages 26 and 31)


Printouts with a caption (page 40)


In addition to the above printouts of multiple reduced images, printouts of multiple reduced images with white borders can be made.

Printouts of identical images (Page 28)

System Configuration

The following shows an example printer system configuration.

Before Printing

This section describes the following operations that must be made prior to start printing after installing the printer and making connections.


- Loading an ink ribbon cassette (see page 8)
- Loading paper (see page 10)
- Selecting the input signal (see page 12)

Once the above operations are done, there should be no need to subsequently perform in routine printing operations. Perform the above operations, if necessary.


Loading an Ink Ribbon Cassette

To make printouts, an ink ribbon cassette and paper should be loaded. Both of those should be used in correct pairs. (see "Ink Ribbon Cassette and Paper" page 90)

- 1** Push in the PUSH OPEN button.
The front panel opens.


- 2** Remove the ink ribbon cassette by pulling down the EJECT lever.


Note

Never put your hand into the ink ribbon cassette dock. The thermal head becomes very hot. You may burn yourself if you touch it.


- 3** Take up any slack in the ink ribbon.
If the ribbon is left slack, it may be crumpled and damaged when inserted.


Wind the spools of the ink ribbon cassette as illustrated until a black bar extending full-width of the ribbon appears on the ink ribbon.

- 4** Insert the ink ribbon cassette firmly until it stops.

- 5** Close the front panel.

Before Printing (continued)

Notes

When using ink ribbon cassettes:


- Once an ink ribbon has been completely used, replace it. Ink ribbons are not reusable.
- Do not touch the ribbon or place the cassette in a dusty place. Body oils or dust stuck to the ink ribbon will cause imperfect printing.

When storing ink ribbon cassette:

- Avoid placing the ink ribbon cassette in a location subject:
 - high temperatures
 - high humidity
 - excessive dust
 - direct sunlight
- Store a partially used ink ribbon in its original bag.

If your ink ribbon should tear

Repair the tear with transparent tape. There should be no problem in using the remaining portion of the ribbon.


Loading Paper

Follow these steps to load paper in the printer. Use only the ink ribbon cassette and paper packed in the same carton, that is correctly in pairs. Be careful not to touch the printing surface.

Note

When loading the paper while the printer is operating, do not turn off the power. If you turn off the power, the image stored in memory will be lost.


- 1 Push PUSH on the paper tray.
The paper tray is ejected.


2 Place the paper into the paper tray.

Notes

- The paper tray holds up to 100 sheets. When you add paper to a partly-full tray, be careful that the total number of sheets does not exceed 100. If you exceed this limit, paper jams may occur.
- Load the paper so that it lays flat in the paper tray.
If the paper is curled, it will overflow the paper tray and the printing position may shift. If this happens, load fewer sheets in the paper tray.

3 Slide the paper tray back into the printer until it clicks into place.

Notes


When storing paper:

- Avoid placing the paper subject to:
 - high temperatures
 - high humidity
 - excessive dust
 - direct sunlight
- Keep the package for storing unused paper.

Before Printing (continued)


Selecting the Input Signal

Before printing, select the input signal (the input connector to which the video signal to be printed is being input (VIDEO, S-VIDEO, or RGB)), from the menu on the video monitor or printer window display. Once you have selected the input signal, this setting remains as is until you select another source.


- 1 Turn on the video monitor and the printer.
The following message appears when the printer is ready to operate.


Video monitor screen


Printer window display

- 2 Press the MENU button.
The following screen appears.


Main Menu screen

Part of the main menu

3 Select INPUT SEL by pressing the ^ or v button.


Move the cursor to INPUT SEL by pressing the ^ or v button.

INPUT SEL :
<VIDEO> S / R

Display this message by pressing the ^ or v button.

4 Select the desired input signal by pressing the < or > button.

Switch the desired input signal to green by pressing the < or > button.
The name of the selected input signal appears in green.

INPUT SEL :
<VIDEO> S / R

Position < > over the name of the desired input signal.

Source signal of the Image to be printed	Video monitor and printer window display (The name of the selected input signal appear on the screen.)
Signal from the video equipment connected to the VIDEO INPUT connector	V → VIDEO
Signal from the video equipment connected to the S-VIDEO INPUT connector	S → S-VIDEO
Signal from the video equipment connected to the RGB/SYNC INPUT connectors	R → RGB


5 Press the MENU button.
The regular screen appears.

Making Full-Size Printouts

This section explains how to make a full-size printout. The operations described here is the basic procedure for making a printout.


Before making a full-size printout

- All connections should have already been made. (see page 52)
- Ensure that the appropriate ink ribbon cassette/paper set is being used and that they are correctly loaded. (see pages 8, 10 and 90)
- Select the input signal to be used to make a printout. (see page 12)
- Set the memory mode to store one full-size image into memory. (see page 27)
- Select the appropriate memory page. (see page 25)


- 1 Turn on the video monitor and the printer.
The following message appears when the printer is ready to operate.


Video monitor screen


Printer window display

- 2** Start the video source.
(This operation is done using the controls of the video equipment acting as the source.)

- 3** Press the MEMORY IN button at the instant when the image you want to print appears on the screen.
That image is stored into memory.
The memory image (stored into memory) is displayed on the screen.

If the stored image is blurred

A quickly moving image may be blurred when it is printed. If this happens, switch the FRM/FLD (frame/field) mode setting to FLD on the main menu and perform printing again. This should eliminate blur from the printout. However, since printing in field mode has a lower resolution than in the frame mode, the ultimate print quality will be slightly degraded. (see "About Memory" page 22)

To change the stored image

- ① Press the SOURCE/MEMORY button.
The image from the video source appears.
- ② Press the MEMORY IN button at the instant when the image you want to print appears.
The previous image is replaced.

Note

If you turn off the power, the image stored into memory will be lost. Thus, store the image into memory again when you turn on the power.

Continue to next page →


Making Full-Size Printouts (continued)

4 Press the PRINT button.

It takes about 60 seconds to make a color printout, or 30 seconds to make a black and white printout.

Blinks while printing.

Does not blink while printing

During color printing: Printing start - yellow - magenta - cyan - printing end
During black and white printing: Printing start - white - printing end

Note

Do not pull the paper from the paper cover until printing has been completed.

To stop printing before completion

Press the STOP button. Printing is abandoned and the paper is ejected to the print tray.

If the printer does not print

When the printer does not print even if the PRINT button is pressed, a message is displayed on the video monitor and printer window display. (see "Error Messages" page 93)

If a black line appears on the printout

Sometimes, a black line appears on the printout, although it does not appear on the video monitor. You can eliminate the black line from the printout. (see "Changing the Printout Area" page 63)

Notes

When preserving your printouts:


- Keep printouts in a dark and cool place.
- Do not stick plastic tape to the print. Also avoid leaving plastic eraser on top of the printout or putting the printout between things which contain plasticizer (a desk mat, etc.).
- Do not pour alcohol or other volatile organic solvents on the printouts.

Queuing images to be printed out

You can store other images to another memory to be printed out once the printer becomes free.

When using the UP-1800EPM printer, select field mode to queue an image for printing.


- 1 Select the memory page to be printed, by pressing the MEMORY PAGE button.

The available memory pages appear in white.

- 2 Press the MEMORY IN button at the instant when the image you want to print appears on the screen.

- 3 Press the PRINT button.
The image selected in step 2 is queued. The image is printed as soon as all previous printing jobs have been completed.

Memory page whose image have been queued for printing (lights in green)
The memory page display returns to white once printing has been completed.

Note

Another image cannot be stored into a memory page whose image has been queued for printing.


- 4 To queue another memory page, repeat steps 1, 2 and 3.
(only for up-1850 EPM)

Making Full-Size Printouts (continued)

Making Multiple Copies of Identical Image


You can print up to 100 copies of a stored image.

Do the following steps before you start printing or while printing. You can change the designated number of copies any time during printing.

- 1** Press the MENU button.
The following screen appears.

Video monitor screen

Main Menu screen


Printer window display

PRINT QTY :

1

Part of the main menu

- 2** Select PRINT QTY by pressing the ^ or v button.

Move the cursor to PRINT QTY by pressing the ^ or v button.


PRINT QTY :

1


Display this message by pressing the ^ or v button.

- 3** Set the number of copies by pressing the < or > button.

When setting	Button
To decrease the quantity	<
To increase the quantity	>

- 4** Press the MENU button.
The regular screen appears.

When paper runs out during printing

Fill the paper tray with paper and press the PRINT button again. (see "Loading Paper" page 10)

Designating the number of copies by the remote control unit (supplied)


You can designate the number of copies directly on the regular screen by using PRINT QTY +/- buttons.

When setting	Button
To decrease the quantity	PRINT QTY -
To increase the quantity	PRINT QTY +

Making Variations of Printouts


You can store various kinds of images into memory by changing the memory mode and can vary the printout of the stored images by changing the print mode. This section explains how to set the memory mode and change the print mode.

Types of images that can be stored into memory

Types of printouts that the printer can produce

By varying the print mode, the following types of printout can be made using images stored in memory.


You can adjust the size of the printouts. (see page 85)

However, the size of a printout of images from multiple memory pages and the size of a printout of multiple reduced image with white borders cannot be changed.

Printing Operation Flowchart

The following flowchart shows the flow of a printing operation.

Making Variations of Printouts (continued)

About Memory

To make a printout, it is first necessary to store the desired image into memory. The method of storing images into memory is called memory mode.

By setting memory mode, you can store a full-size image or multiple reduced images into memory. Also, the printer can automatically store images to its memory pages at preset intervals.

Also, you have to decide how to use the printer's memory to store images. Two methods of using memory are supported. One is frame mode, while the other is field mode. The number of memory images you can store depends on whether you select frame or field mode.

Frame mode and field mode

Frame (FRM) mode

A image for one screen is stored in one memory.

Video monitor

to 1/3 frame memory page

Printout

1/3FRM: The first page of three frame-mode memory pages

A still subject can be printed with high resolution.


Field (FLD) mode

One memory is divided into two, and images for the two screens are stored to the resulting memory pages.

Video monitor

to 1/6 field memory page

Printout

1/6FLD: The first memory page of six field-mode memory pages

A quickly moving subject can be printed without blurring.

The above illustrations show the case of UP-1850EPM. In case of UP-1800EPM, there is one memory page in frame mode and two memory pages in field mode.

Making Variations of Printouts (continued)

About Memory

To make a printout, it is first necessary to store the desired image into memory. The method of storing images into memory is called memory mode.

By setting memory mode, you can store a full-size image or multiple reduced images into memory. Also, the printer can automatically store images to its memory pages at preset intervals.


Also, you have to decide how to use the printer's memory to store images. Two methods of using memory are supported. One is frame mode, while the other is field mode. The number of memory images you can store depends on whether you select frame or field mode.

Frame mode and filed mode


Frame (FRM) mode

A image for one screen is stored in one memory.


Video monitor

to 1/3 frame memory page

Printout


1/3FRM: The first page of three frame-mode memory pages

A still subject can be printed with high resolution.


Field (FLD) mode

One memory is divided into two, and images for the two screens are stored to the resulting memory pages.

Video monitor

to 1/6 field memory page

Printout


1/6FLD: The first memory page of six field-mode memory pages

A quickly moving subject can be printed without blurring.

The above illustrations show the case of UP-1850EPM. In case of UP-1800EPM, there is one memory page in frame mode and two memory pages in field mode.


Selecting frame or field mode

Before storing an image, select frame or field mode.


- 1** Press the MENU button.
The following screen appears.

Video monitor screen


Main Menu screen

Printer window display

Part of the main menu

- 2** Select FRM/FLD by pressing the \wedge or \vee button.

Move the cursor to FRM/FLD by pressing the \wedge or \vee button.


Press the \wedge or \vee buttons until FRM/FLD appears.

Continue to next page →

Making Variations of Printouts (continued)

- 3** Select the desired mode by pressing the < or > button.

Switch the desired mode to green by pressing the < or > button.

Position < > over the desired mode by pressing the < or > button.

FRAME: We recommend that, whenever possible, you print in this mode.

FIELD: Select this mode to reduce blurring when you print a quickly moving image.

- 4** Press the MENU button.
The regular screen appears.

About memory pages

The UP-1800EPM has a single frame memory, enabling the unit to store one image in one memory page when FRM mode is selected, or two images in two memory pages when FLD mode is selected.



The UP-1850EPM has three frame memories, enabling the unit to store three images in three memory pages in frame mode and six images in six memory pages in field mode.

The memory used to store one screen image is called a memory page.

Selected memory mode	Number of usable memory pages		Usable memory pages	
	UP-1800EPM	UP-1850EPM	UP-1800EPM	UP-1850EPM
Frame mode (FRM)	1	3	1/1FRM	1/3FRM, 2/3FRM or 3/3FRM
Field mode (FLD)	2	6	1/2FLD 2/2FLD	1/6FLD, 2/6FLD, 3/6FLD, 4/6FLD, 5/6FLD or 6/6FLD

Selecting a memory page

To select a memory page, press the MEMORY PAGE button.

Press the MEMORY PAGE button as many times as necessary until the desired memory page appears.

Making Variations of Printouts (continued)

Selecting the Memory Mode


Decide the method for storing images in memory. Once you have selected memory mode, this setting remains as is until reset, even if you turn the power off.

This section explains the following settings.

- Selecting the number of images stored in one memory page (full-size image/multiple reduced images)
- Automatic storing of images into a memory page at preset intervals

When you use the remote commander to control the printer

You can access the MULTI PICTURE sub menu by pressing the MULTI PICTURE button. Thus, press the MULTI PICTURE button to display the MULTI PICTURE sub menu. Then, follow the procedure below, starting from step 4.

- 1 Press the MENU button.
The following screen appears.

Video monitor screen

```
MENU                                COLOR 100
▶ PRINT QTY : 1
  INPUT SEL : VIDEO/S/R
  FRM/FLD   : FRAME/FIELD
  SHARPNESS : L/M/H
  MULTI PIX : PUSH [ > ]
  COLOR ADJ : PUSH [ > ]
  CAPTION   : PUSH [ > ]
  SET UP    : PUSH [ > ]
```

Main Menu screen

Printer window display

```
PRINT QTY : 1
```

Part of the main menu

- 2 Select MULTI PIX by pressing the ^ or v button.


```
MENU                                COLOR 100
  PRINT QTY : 1
  INPUT SEL : VIDEO/S/R
  FRM/FLD   : FRAME/FIELD
  SHARPNESS : L/M/H
  MULTI PIX : PUSH [ > ]
  COLOR     : PUSH [ > ]
```

Move the cursor to MULTI PIX by pressing the ^ or v button.


```
MULTI PIX :
          PUSH [ > ]
```

Press the ^ or v buttons until MULTI PIX appears.

- 3** Press the > button.
The following screen appears.


MULTI PICTURE sub menu



Part of the MULTI PICTURE sub menu

- 4** Set the memory mode.
- ① Select the item to be set by pressing the ^ or v button.
 - ② Select the method for storing images by pressing the < or > button.

Switch the desired mode to green by pressing the < or > button.

Position <> over the desired mode by pressing the < or > button.

When you select	Item for memory mode	Settings	Contents of setting
To set the number of images to be stored in one memory page.	SPLIT MEM	OFF	Storing a full-size image
		4	Storing four reduced images
		16	Storing 16 reduced images
To set whether the printer automatically stores images into memory page at preset intervals.	INTERVAL	OFF	Storing only the image currently displayed on the screen by pressing the MEMORY IN button.
		1/2/3/4	Storing images sequentially to all memory pages at preset intervals. (1: about 1/60 seconds, 2: about 1/15 seconds, 3: about 1/4 seconds and 4: about 1/2 seconds)

- 5** Press the MENU button.
The regular screen appears.

Continue to next page →

Making Variations of Printouts (continued)


Selecting the Print Mode

You can make variations of printouts from the images stored in memory pages by changing the print mode. (see “Types of printouts that the printer can produce” page 20)

This section explains the following settings.


- Deciding whether images are to be printed with white borders (printout with/without white borders)
- How many identical images are printed on a single printout
- Setting the size of the image in a single printout
- Printing an image rotated around its vertical axis (mirror image printout)

Once you have selected the print mode, this setting remains as is until you reset, even if you turn the power off.


- 1** Press the MENU button.
The following screen appears.

Video monitor screen


Main Menu screen

Printer window display

Part of the main menu

2 Select MULTI PIX by pressing the ^ or v button.


Move the cursor to MULTI PIX by pressing the ^ or v button.

MULTI PIX :
PUSH [>]

Press the v or ^ buttons until MULTI PIX appears.

3 Press the > button. The following screen appears.


MULTI PICTURE sub menu

SPLIT MEM :
<OFF> 4 / 16

Part of the MULTI PICTURE sub menu

4 Set the print mode.

- ① Select the item to be set by pressing the ^ or v button.
- ② Select the method for making a printout by pressing the < or > button.

Switch the desired mode to green by pressing the < or > button.

SEPARATE :
<OFF> ON

DUP IMAGE :
<OFF> 2 / 4 / 16

IMAGE SIZE :
<L> M / S

MIRROR :
<OFF> ON

Position <> over the desired mode by pressing the < or > button.

Some items affect the setting of other items. That is, there are some items whose setting displays will change when the setting of another item is changed.

Continue to next page →

Making Variations of Printouts (continued)

When you select	Item for print mode	Settings	Content of settings
To decide whether the images are printed with white borders	SEPARATE	OFF	without white borders ^{a)}
		ON	with white borders ^{b)}
To decide how many times identical reduced images are printed in a single printout.	DUP IMAGE ^{a)}	OFF	Printing a memory image one time
		2	Printing a memory image twice
		4	Printing a memory image four times
		16	Printing a memory image 16 times.
To set the size of the image in a single printout.	IMAGE SIZE ^{a)}	L	Regular size
		M	About 70% of the L size
		S	About 50% of the L size
To rotate the image around its vertical axis (to make a mirror image printout)	MIRROR	OFF	Normal image
		ON	Mirror image
To make a printout consisting of multiple images, each from a same memory page or multiple memory pages.	PICK IMAGE	This item is effective only when the PICK sub menu is displayed.	

a) When SEPARATE is set to OFF (without white borders), the following items can be set.

DUP IMAGE	IMAGE SIZE		
	L	M	S
OFF	O	O	X
2	O	O	X
4	O	O	O
16	O	X	X

O: Size that can be set
X: Size that cannot be set

b) White borders can be added only for the following settings.

SPLIT MEM: Settings other than OFF

DUP IMAGE: OFF

IMAGE SIZE: L

- 5 Press the MENU button.
The regular screen appears.

To return to the main menu from the sub menu screen


Press the EXEC button except when the SAVE item is selected on the COLOR ADJUST sub menu and when the cursor is position in the character entry area on the CAPTION sub menu.

Making Printouts of Multiple Reduced Images


This subsection explains how to make printouts of multiple reduced images taking as an example, making a printout of 16 reduced images. (see "Selecting the Print Mode" page 28)

Before making printouts of 16 reduced images

- Set the memory mode to store 16 reduced images into memory. (see page 27)
- Select the appropriate memory page. (see page 25)

- 1 Start the video source.
(This operation is done using the controls of the video equipment acting as the source.)

Continue to next page →


Making Variations of Printouts (continued)

- 2** Press the MEMORY IN button at the instant when the image you want to print appears on the screen.

The image is stored to the position for which the corresponding number blinks on the monitor display, and over which the corresponding number < > is positioned in the printer window display.

The cursor moves to the next number.

When an image has been stored into this position, the number lights in green.


- 3** Press the SOURCE/MEMORY button.
The image from the video equipment appears on the monitor display.
- 4** Repeat steps 2 and 3 until you have stored 16 images.

To change a stored image

Example: When you want to change the image stored to the 5th position.

- ① Select 5 by pressing the ^, v, < or > button.

- ② Press the SOURCE/MEMORY button.
The image from the video source appears.
- ③ Press the MEMORY IN button at the instant when the image you want to print appears.
The previously stored image is replaced with the newly selected image.

To skip a previously stored image

When an image has already been stored, the previously stored image can be replaced by pressing the MEMORY IN button.

Skip the number corresponding to the image to be skipped by pressing the ^, v, < or > button.

5 Set the print mode. (see “Selecting the Print Mode” page 28)

6 Press the PRINT button.
The 16 reduced images are printed on one sheet of paper.

To stop printing midway

Press the STOP button. The printer stops printing and ejects paper to the paper cover.

To automatically store images into memory

Images can be automatically stored into all memory pages at preset intervals by using the INTERVAL item of the MULTI PICTURE sub menu, which is used to set the memory mode. (see “Selecting the Memory Mode” page 26)

While images are automatically being stored to all memory pages, “PLEASE WAIT” appears on the screen.

Making printouts with an insert

You can make printouts with an insert by using the four- or 16- reduced image memory mode.

Example: To make a printout with one of four reduced images inserted



- 1** Display the full-size image stored in memory. (Follow steps 1 to 3 of “Making Full-Size Printouts” on page 14.)
- 2** Set the memory mode to store four reduced images. (see “Selecting the Memory Mode” page 26)

Continue to next page →

Making Variations of Printouts (continued)

- 3** Move the white blinking cursor to the position where a reduced image is to be inserted, by pressing the ^, v, < or > button.
Example: To insert the image to 2

Move the white blinking cursor to 2.

- 4** Press the SOURCE/MEMORY button to display the image from the video source, if necessary.
- 5** Press the MEMORY IN button at the instant when the image you want to print appears.
The image is stored to position 2.
- 6** Press the PRINT button.
An image with the insert is printed.


Note

If you insert a reduced image into an image stored in a different memory page, the printer can not make a printout of the image with an insert.

Making Printouts of Multiple Images in Multiple Memory Pages

You can make a single printout of multiple images stored in multiple memory pages. The number of images of a single printout depends on the memory mode setting (the setting of SPLIT MEM from the MULTI PICTURE sub menu). In this manual, the type of printout described above is referred to as a composed image.

Type of memory images	Number of images in one composed image printout
Full-size image	2
Four reduced images	4
16 reduced images	16

When using the UP-1800EPM, which has only a single frame memory, select field mode to make a printout of multiple images in multiple memory pages.

You can also exchange the position of multiple reduced images in the same memory page.


Before making printouts of images in multiple memory pages (composed image)

- Store images of the same type to each memory page.
- Display the images stored in memory on the screen.

Note


Store the images of the same type to each memory page (one full-size image, four reduced images or 16 reduced images). Otherwise, the printout is likely to be unsatisfactory. For example, there may be whitish areas or different images may be overlapped.

Making Variations of Printouts (continued)


- 1** Press the MENU button.
The following screen appears.

Video monitor screen


Main Menu screen

Printer window display

Part of the main menu

- 2** Select MULTI PIX by pressing the \wedge or \vee button.


Move the cursor to MULTI PIX by pressing the \wedge or \vee button.

```

MULTI PIX :
PUSH [>]
  
```

Press the \wedge or \vee buttons until MULTI PIX appears.

- 3** Press the > button.
The following screen appears.



MULTI PICTURE sub menu

```

SPLIT MEM :
<OFF> 4 / 16
  
```

Part of the MULTI PICTURE sub menu

- 4** Select the type of the composed image.
- ① Select SPLIT MEM by pressing the \wedge or \vee button.
 - ② Select the number of images in the composed image by pressing the $<$ or $>$ button.


Position $<>$ over the desired number.


Switch the desired number to green.

Selection from SPLIT MEM	Number of images in one composed image printout
OFF	2
4	4
16	16

- 5** Select PICK IMAGE by pressing the \wedge or \vee button.


Move the cursor to PICK IMAGE by pressing the \wedge or \vee button.

Press the \wedge or \vee button until PICK IMAGE appears.

- 6** Press the $>$ button.
The following screen appears.

PICK sub menu

Part of the PICK sub menu

Blinks: shows that the image will be inserted into A.

To store an image to a position other than A


Images are stored automatically in order of A, B, C, D..... You can select where each image is to be stored.

- ① Select PRINT POSITION by pressing the \wedge or \vee button.
- ② Select the desired position by pressing the $<$ or $>$ button.

Continue to next page →

Making Variations of Printouts (continued)

- 7** Select MEMORY POSITION by pressing the \wedge or \vee button.


Move the cursor to POSITION by pressing the \wedge or \vee button.

MEM POS : 1..4
<1> 2

Press the \wedge or \vee button until MEM POS appears.


- 8** Select the memory page to which the desired image is stored, by pressing the MEMORY PAGE button.

The memory page in which the image displayed on the video monitor is stored.

MEM POS : 1..4
<1> 2


- 9** Select the image to be stored by pressing the < or > button.

MEM POS : 1..4
<1> 2

The selected image (that to which the cursor is positioned) blinks. This image will be stored to the position that is blinking in the PRINT POSITION item of the PICK sub menu.

- 10** Press the EXEC button.
The image selected in step 9 is stored.

MEM POS : 1..4
<1> 2

Blinks (in the next step, the image to be stored here will be decided).

Lights in green (the image is stored).

11 Repeat steps 8, 9 and 10 until you have stored all the desired images.

To leave a position blank

A position can be left blank without inserting an image.

- ① Select PRINT POSITION by pressing the ^ or v button.
- ② Select the position that you want to leave blank, by pressing the < or > button.
- ③ Press the EXEC button.
Any image is not printed at that position.

To change a stored image

- ① Select PRINT POSITION by pressing the ^ or v button.
- ② Select the position of the image to be changed by pressing the < or > button.
- ③ Select MEMORY POSITION by pressing the ^ or v button.
- ④ Store a new image.
The previous image is replaced with the new one.

To check an image that is hidden below a screen message

You can check an image by temporarily deleting screen display.

- ① Select CHECK by pressing the ^ or v button.
- ② Press the > button.
For as long as you keep the > button held down, the display does not appear on the screen.

12 Press the PRINT button.

The composed images are printed on a single paper.

Note

To make a printout of a composed image, press the PRINT button on the PICK sub menu. The composed image can not be printed on the regular screen or main menu screen.

Making Printouts with a Caption

A caption, such as data or comments, can be added to a printout, using small characters below the image.

You can input up to 128 characters in NARROW size mode or NORMAL size mode, and up to 64 characters in WIDE size mode.

Note

When the printout is printed in field mode, characters may not be printed clearly.

About the CAPTION sub menu


A caption is entered from the CAPTION sub menu. A brief explanation of each item on the CAPTION sub menu, is given below before entering a caption.

Character display area

CAPTION ON: displayed when printing with a caption
CAPTION OFF: displayed when printing without a caption
MIRROR: displayed when printing with mirror characters

The position where the cursor is placed in character display area (the position where the character is to be entered)

The line where the cursor is positioned in the character display area is displayed.

SAMPLE 1 2 3
BS SP OF ON MI SAVE

The line where the cursor is positioned in the character entry area is displayed.

Symbols and words with = are used to enter a caption.

Character entry area

The cursor is positioned at the highlighted character and this highlighted character is to be entered.

Note

Lower-case ä, ö and ü are displayed instead of capital Ä, Ö and Ü even in the capital-letter setting and ≈ is displayed instead of ~ on the printer window display. However, in actual printouts, those letters and the symbol are printed correctly according to the settings.


Symbols and words with = used to enter a caption

Monitor display	Printer window display	Function
SPACE	SP	One space
BS	BS	One backspace
OFF	OF	Selecting to print without a caption
ON	ON	Selecting to print with a caption
MIR	MI	Selecting to print with a mirror caption
SHIFT ^{a)}	SF	Selecting either capital letters or lower-case letters
SAVE	SAVE	Storing the entered caption

a) By highlighting SHIFT and pressing the EXEC button, capital letters are changed to lower-case letters, or lower-case letters are changed to capital letters.


Entering a Caption

Enter a caption as follows. The setting remains valid until you enter a new setting - even if you turn the power off.


- 1 Press the MENU button.
The following screen appears.

Video monitor screen

Main menu screen

Printer window display


Part of the main menu

Continue to next page →

Making Printouts with a Caption (continued)

- 2** Select CAPTION by pressing the \wedge or \vee button.


Move the cursor to CAPTION by pressing the \wedge or \vee button.


CAPTION:
PUSH [>]

Press the \wedge or \vee buttons until CAPTION appears.

- 3** Press the > button.
The following screen appears.



CAPTION sub menu


Part of the CAPTION sub menu

- 4** Select the position where you want to enter the character in the character display area by pressing the < or > button.


The cursor \square is highlighted at the selected position on the monitor display.
The cursor \square blinks in black at the selected position on the printer window display.

- 5** Select the character you want to enter by pressing the \wedge , \vee , < or > button.
Example: To select S

Highlight S.

Press the \wedge or \vee button until this line appears, then press < or > buttons until S blinks.

6 Press the EXEC button.

The selected character appears at the position highlighted on the character display area, then the highlighted  moves to the next position.

When you enter a wrong character

Select BS by pressing the ^, v, < or > buttons, then press the EXEC button. The character to the left of highlighted character will be deleted.

7 Repeat steps 4, 5 and 6 to enter the remaining characters of the caption.

To make a space

- ① Move the highlighted  to the position where you want to make a space.
- ② Select SPACE by pressing the ^, v, < or > button.
- ③ Press the EXEC button.

The one space is made and the cursor moves to the next position.


To replace a previously entered character without changing the number of characters

You can replace a previously entered character with a new one.

- ① Move the cursor to the character which you want to replace by the operation in step 4.
- ② Enter the correct character over the wrong character by the operations in steps 5 and 6.

The previously entered character is replaced with the new one.

8 Select SAVE by pressing the ^, v, < or > button.

Press the ^ or v button until this line appears, then press the < or > button until SAVE blinks.

9 Press the EXEC button.

The message "PLEASE WAIT" appears while the entered characters are being stored. Once they have been stored, the message disappears and the CAPTION sub menu appears again.



10 Press the MENU button.
The regular screen appears.

Making Printouts with a Caption (continued)

Making printouts with a caption

Display the CAPTION input screen. (see “Entering a Caption” page 41)

- 1 Select ON by pressing the ^, v, < or > button.

Highlight ON.

- 2 Press the EXEC button.

Making a printouts without a caption

Select OFF in the above step 1.

Making a printout with a mirror caption


You can also rotate a caption around its vertical axis (mirror caption) when making a mirror image printout.

Check

The MIRROR setting of the MULTI PICTURE sub menu should be set to MIRROR ON. (see page 28)

Display the CAPTION input screen. (see “Entering a Caption” on page 41)

- 1 Select MIR by pressing the ^, v, < or > button.



Highlight MIR.

When the error tone sounds three times

The MIRROR setting of the MULTI PICTURE sub menu is set to MIRROR OFF.

- 2 Press the EXEC button.

MIRROR is displayed.

To return the print mode to the one with normal caption

- ① Select ON by pressing the ^, v, < or > button.
- ② Press the EXEC button.

To return to the regular screen

Press the MENU button.

Deleting the Images Stored into Memory Pages


You can delete images stored to memory pages, either individually or all at once, from a menu.

Note

You cannot restore images once they have been deleted.


Simultaneously Deleting All Images in a Single Memory Page

You can delete a full-size image or four or 16 reduced images from one memory page.


- 1** Press the SOURCE/MEMORY button when the image from the video equipment appears on the monitor display.
- 2** Press the MENU button.
The following screen appears.

Video monitor screen

Main menu screen

Printer window display


Part of the main menu

Continue to next page →

Deleting the Images Stored into Memory Pages (continued)

- 3** Select SET UP by pressing the \wedge or \vee button.


Move the cursor to SET UP by pressing the \wedge or \vee button.

SET UP :
PUSH [>]

Press the \wedge or \vee button until SET UP appears.

- 4** Press the > button.
The following screen appears.


SET UP sub menu

H SHIFT :
0 DOT

Part of the SET UP sub menu

- 5** Select MEM CLR by pressing the \wedge or \vee button.


Move the cursor to MEM CLR by pressing the \wedge or \vee button.

MEM CLR :
PUSH [>]

Press the \wedge or \vee button until MEM CLR appears.

- 6** Press the > button.
The following screen appears.

MEMORY CLEAR sub menu


ALL CLR :
PUSH [EXEC]

Part of MEMORY CLEAR sub menu


When you want to see an image that is hidden behind the menu
You can erase the screen message.

- ① Select CHECK by pressing the \wedge or \vee button.
- ② Press the > button.
For as long as you keep the > button held down, the menu does not appear on the screen.

- 7** Select ALL CLR by pressing the \wedge or \vee button.

Move the cursor to ALL CLR by pressing the \wedge or \vee button.


Press the \wedge or \vee button until ALL CLR appears.

- 8** Select the memory page, in which the images to be deleted are stored, by pressing the MEMORY PAGE button.
- 9** Press the EXEC button.
All images in the selected memory page are deleted.
- 10** Repeat steps 8 and 9 to delete the images in the remaining memory pages.
- 11** Press the MENU button.
The regular screen appears.

Deleting the Images Stored into Memory Pages (continued)

Deleting Part of Multiple Images in a Single Memory Page


You can delete unnecessary images from a single memory page.

- 1** Press the SOURCE/MEMORY button when the source image from the video equipment appears on the monitor.
The memory image appears on the monitor.

- 2** Press the MENU button.
The following screen appears.

Video monitor screen


Main Menu screen

Printer window display


Part of the main menu

- 3** Select SET UP by pressing the ^ or v button.


Move the cursor to SET UP by pressing the ^ or v button.


Press the ^ or v button until SET UP appears.

- 4** Press the > button.
The following screen appears.


SET UP sub menu


Part of the SET UP sub menu

- 5** Select MEM CLR by pressing the ^ or v button.


Move the cursor to MEM CLR by pressing the ^ or v button.


Press the ^ or v button until MEM CLR appears.

- 6** Press the > button.
The following screen appears.

MEMORY CLEAR sub menu

Part of MEMORY CLEAR sub menu


When you want to see an image that is hidden behind the menu
You can erase the screen message.

- ① Select CHECK by pressing the ^ or v button.
- ② Press the > button.
For as long as you keep the > button held down, the menu does not appear on the screen.

Continue to next page →

Deleting the Images Stored into Memory Pages (continued)

- 7** Select PART CLR by pressing the \wedge or \vee button.


Move the cursor to PART CLR by pressing the \wedge or \vee button.

PART CLR : 1..16
<1> PUSH [EXEC]

Press the \wedge or \vee button until PART CLR appears.

- 8** Select the memory page containing the images to be deleted stored by pressing the MEMORY PAGE button.

- 9** Select the image to be deleted by pressing the < or > button.
Example: Deleting the sixth image of 16 reduced images

PART CLR : 1..16
<6> PUSH [EXEC]

Press the < or > button until 6 appears.

- 10** Press the EXEC button.
The image selected in step 9 is deleted.


- 11** Repeat steps 9 and 10 to delete the remaining images.

- 12** To delete images in other memory pages, repeat steps 8, 9, 10 and 11.


- 13** Press the MENU button.
The regular screen appears.

Supplied Accessories


Remote commander RM-5100 (1)


Paper tray (1)


Connecting cable for remote commander (1)


Paper cover (1)

Dry batteries (2)

AC power cord (1)

Color printing pack UPC-1010 (1)
Instructions For Use (1)

Connections

To enable printing, video equipment to act as an input signal source, and a video monitor to enable you to view images or menus, must be connected. The following diagrams illustrate how to make the input, output and remote control connections. Use as a guide when connecting the necessary signals to and from the equipment to be used for printing.


Notes

When connecting:

- Turn off the power of each device before attempting to make any connections.
- Connect the AC power cord last.

Making Connections for Storing Video Images

Connect the video equipment for storing the video images to be printed. Connect the necessary video equipment which will be used in actual printing, using the following diagram as a guide. Before connecting the video equipment, see "Important safeguards/notices for use in the medical environment" on page 2.


a) Normally, set this switch to ON. Set it to OFF if the level of the input signal drops if you connect additional video equipment.

Making Connections for Viewing Images to be Printed on the Video Monitor

Connect a video monitor to view stored images and to check those to be printed. Connect the necessary video monitor which will be used in actual printing, using the following diagram as a guide.

Before connecting the video monitor, see "Important safeguards/notices for use in the medical environment" on page 2.

a) When you connect the video monitor to RGB/SYNC OUTPUT, set SYNC ON G of the MONITOR ADJUST sub menu to ON.

To set SYNC ON G to ON

- 1** To access the MONITOR ADJUST sub menu, perform steps 1 to 5 of "Matching the Video Monitor Color to the Printer Color" on page 69.
- 2** Select SYNC ON G by pressing the \wedge or \vee button.
- 3** Set SYNC ON G to ON by pressing the $<$ or $>$ button.
- 4** Press the MENU button.
The regular screen appears.

Connections (continued)


Making Connections to Enable Remote Control

The printer can be controlled remotely by connecting the remote control unit (supplied), the RM-91 remote commander (see "Preparing the Remote Control Units" page 55), or a personal computer.

Before connecting the computer, see "Important safeguards/notices for use in the medical environment" on page 2.

Note


Connect the supplied remote commander to REMOTE 1 on the front panel. If you accidentally connect the supplied remote commander to REMOTE 2, it may be difficult to disconnect.

- a) When connecting a personal computer, select the appropriate baud rate by setting the DIP switches SW3 and SW4. (see "Setting the DIP switch" page 85)

Assembly

Mount the supplied paper tray and paper cover.

Preparing the Remote Commanders

You can control the printer remotely by using the remote commander (supplied) or the remote commander (not supplied).


Using the Supplied Remote Commander

The remote commander can be used either as a wireless type or wired type. The buttons on the remote commander duplicate those on the front panel of the printer, except for the PRINT QTY button, COLOR ADJUST button and MULTI PICTURE button. (see "Remote Commander RM-5100" page 101)


Inserting batteries

Install the batteries in the remote commander before using it as a wireless unit.

- 1 Remove the battery compartment cover.

- 2 Insert the two supplied SUM-3 1.5 V batteries. Note the polarity. Be careful to insert the batteries correctly.

- 3 Replace the cover.

Battery life

The battery life depends on how much you use the remote commander. On average, batteries last for about 6 months. Install fresh batteries as soon as you notice the unit's range becoming shorter.

Notes


When using the batteries:

- Remove the batteries from the remote commander if you do not intend to use it for an extended period of time. The batteries may leak if you leave them in the remote control unit.
- Should the batteries leak, clean the battery case thoroughly with a soft cloth and install fresh batteries.
- Be careful to insert the batteries correctly. Note the polarity, as indicated inside the battery compartment.
- Replace exhausted batteries with fresh ones. Never mix a fresh battery with a used battery or with a different kind of battery.

Preparing the Remote Commanders (continued)

Using the supplied remote commander as a wireless unit

When using the remote commander as a wireless unit, aim the head of the remote control unit of the remote sensor on the printer. With fresh batteries, the range of the remote control unit is about 7 meters.

Using the Remote Commander (Not Supplied)

The RM-91 remote commander (not supplied) allows you to make prints remotely.

Operation


At the instant the image you want to print is displayed on the monitor, press the switch of the remote commander. The subsequent operation of the printer will depend on the remote operation setting with the corresponding menu. (see "Selecting the Operation Mode for Automatic Printing Capabilities" page 78) The printer operation, also, can be controlled remotely by sending a pulse signal to the REMOTE 2 connector. (see "Specifications" page 92)

Adjusting the Printout Quality

You can adjust the printout quality, including its sharpness and color (intensity and contrast) and store these settings by using the menu. The setting remains as is until reset - even if you turn off the power.

Adjusting the Sharpness

You can set the printout sharpness to one of three levels: L (Low), M (Medium) or H (High). A printout will appear softer or sharper depending on the definition of the subject outline. The image on the monitor is not affected by changing the sharpness setting. This adjustment affect only the quality of the printout. The setting remains as is until reset - even if you turn off the power.

- 1 Press the MENU button.
The following screen appears.

Video monitor screen

```
MENU                                COLOR 100
▶PRINT QTY : 1
INPUT SEL : VIDEO/S/R
FRM/FLD   : FRAME/FIELD
SHARPNESS : L/M/H
MULTI PIX : PUSH [ > ]
COLOR ADJ : PUSH [ > ]
CAPTION   : PUSH [ > ]
SET UP    : PUSH [ > ]
```

Printer window display

PRINT QTY :

1

Adjusting the Printout Quality (continued)

- 2** Select SHARPNESS by pressing the \wedge or \vee button.

```
MENU          COLOR 100
PRINT QTY    : 1
INPUT SEL    : VIDEO/S/R
FRM/FLD      : FRAME/FIELD
SHARPNESS    : L/M/H
MULTI PIX    : PUSH [>]
```

Move the cursor to SHARPNESS by pressing the \wedge or \vee button.

SHARPNESS :
<L> M / H

Display SHARPNESS by pressing the \wedge or \vee button.

- 3** Select desired sharpness by pressing the < or > button.

```
MENU          COLOR 100
PRINT QTY    : 1
INPUT SEL    : VIDEO/S/R
FRM/FLD      : FRAME/FIELD
SHARPNESS    : L/M/H
MULTI PIX    : PUSH [>]
```

Switch the desired sharpness to green by pressing the < or > button.

SHARPNESS :
<L> M / H

Position <> over the desired sharpness by pressing the < or > button.

Desired sharpness	Content of settings
L (Low)	Soft outline
M (Medium)	Normal outline
H (High)	sharp outline

- 4** Press the MENU button.
The regular screen appears.

Adjusting the Printout Color

This subsection explains how to adjust the printout color.

You can adjust the color intensity (RED/GREEN/BLUE) and contrast (DARK/LIGHT). The new setting remains as is until reset - even if you turn off the power.

You can store up to three settings. These settings are managed according to a LOAD number. The color intensity and picture contrast of a printout are determined by recalling one of the three settings according to their LOAD numbers. The printer retains these settings even if you turn off the power. This is useful when you are using more than one video equipment, each of a different quality, and when you want to print images having different color qualities and picture contrasts.


Also, you can make a printout using temporarily set values, without erasing the stored adjustment values.

Perform the adjustments while viewing the images stored in memory.

All values are factory-set to 0 for LOAD numbers 1, 2 and 3.

When you control the printer using the remote control unit (supplied)


You can directly access the COLOR ADJUST sub menu from the regular screen by pressing the COLOR ADJUST button. Therefore, press the COLOR ADJUST button first. Then, perform the operation from step 5 of the following procedure.

1 Display the image stored in monitor for adjustment.


2 Press the MENU button.
The following screen appears.

Video monitor screen


Main Menu screen

Printer window display


Part of the main menu

3 Select COLOR ADJ by pressing the ^ or v button.

Move the cursor to COLOR ADJ by pressing the ^ or v button.

Press the ^ or v buttons until COLOR ADJ appears.

Continue to next page →

Adjusting the Printout Quality (continued)

- 4** Press the > button.
The following screen appears.

```

COLOR ADJUST
▶LOAD      : 1/2/3
RED        : 0.....
GREEN      : 0.....
BLUE       : 0.....
DARK       : 0.....
LIGHT      : 0.....
CHECK      : PUSH [>]
SAVE       : 1/2/3 [EXEC]
    
```

COLOR ADJUST sub menu

```

LOAD      :
<1> 2 / 3
    
```

Part of the COLOR ADJUST sub menu

- 5** Select LOAD by pressing the ^ or v button.

```

COLOR ADJUST
LOAD      : 1/2/3
RED        : 0.....
GREEN      : 0.....
BLUE       : 0.....
DARK       : 0.....
    
```

Move the cursor to LOAD by pressing the ^ or v button.

```

LOAD      :
<1> 2 / 3
    
```

Press the ^ or v button until LOAD appears.

- 6** Select the LOAD number of the value to be adjusted or to be modified by pressing the < or > button.

```

COLOR ADJUST
▶LOAD      : 1/2/3
RED        : 0.....
GREEN      : 0.....
BLUE       : 0.....
DARK       : 0.....
    
```

Switch the desired LOAD number to green by pressing the < or > button.

```

LOAD      :
<1> 2 / 3
    
```

Position < > over the desired LOAD number by pressing the < or > button.

When modifying, you can preserve the original settings. (see "To preserve the original set value" page 62)

7 Adjust the printout color.

- ① Select the item to be set by pressing the \wedge or \vee button.
- ② Perform the adjustment by pressing the $<$ or $>$ button.


Adjustment item	Contents of setting	
Color intensity	RED	Adjusting the red component of the image
	GREEN	Adjusting the green component of the image
	BLUE	Adjusting the blue component of the image
Color contrast	DARK	Adjusting the dark area of an image
	LIGHT	Adjusting the light area of an image

The RED, GREEN and BLUE color components and the contrast are divided into 15 scales from -7 to $+7$, as indicated by a value and graph. And the center of the graph corresponds to the standard color.

Note


There is an unusable extra scale at the left end on the monitor display.

When adjusting RED/GREEN/BLUE

The intensity increases in the $+$ direction by pressing the $>$ button.
The intensity decreases in the $-$ direction by pressing the $<$ button.

When adjusting DARK/LIGHT

The contrast in the dark area or light area is strengthened in the $+$ direction by pressing the $>$ button. The contrast in the dark area or light area is weakened in the $-$ direction by pressing the $<$ button.

Once you have changed the value

Once you have changed the value, TEMP (TEMPORARY) appears to the right of the LOAD item. TEMP indicates that the setting is temporary and not stored.

Continue to next page →

Adjusting the Printout Quality (continued)

8 After you have made all necessary adjustments, check your presets.


① Select CHECK by pressing the \wedge or \vee button.

② Press the $>$ button.


For as long as you keep the $>$ button held down, the display does not appear on the screen.

You can make a printout with the settings made as above. Go to step 12 to make a printout. However, this setting is cleared when you turn the printer off or you select another preset. To store a new setting, go to the next step.

9 Select SAVE by pressing the \wedge or \vee button.


Move the cursor to SAVE by pressing the \wedge or \vee button.


Press the \wedge or \vee button until SAVE appears.

10 Select the SAVE number to which new settings are to be stored by pressing the $<$ or $>$ button.

Switch the desired SAVE number to green by pressing the $<$ or $>$ button.

Position $<>$ over the desired SAVE number by pressing the $<$ or $>$ button.

To preserve the original set value

Select the SAVE number which is different from the LOAD number selected in step 6.

11 Press the EXEC button.

The settings have been registered to the SAVE number selected in step 10.

When the LOAD number selected in step 6 and the SAVE number selected in step 10 are the same, TEMP disappears from the LOAD item.

12 Press the MENU button.

The regular screen appears.

To recall settings

You can recall previously set values by selecting the LOAD number. The values are stored as SAVE numbers in steps 10 and 11. This SAVE number is the LOAD number for this setting.

Printer Initial Setup

You can set up the following, using the on-screen menu and the DIP switches on the rear panel.

Setup with the menu

- Changing the printout area (see page 63)
- Compensating for input signals (see page 66)
- Matching the video monitor color to the printer color (see page 69)
- Erasing the screen display (see page 75)
- Selecting the operation mode for automatic printing capabilities (see page 78)

DIP switch settings

- Setting the printout size (see page 85)
- Setting the baud rate for computer communication (see page 85)
- MENU button function ON/OFF (see page 85)
- Selecting whether to lock to an external synchronous signal (GENLOCK) (see page 86)
- Selecting whether to enable the operation and error tones (see page 86)

Changing the Printout Area

The black line may be printed on the printout although it does not appear on the video monitor. The portion where no video signal exists is printed in black. This may occur when you make printouts after you connect a different video source or play back different video software.

In such a case, you can adjust the printout area by moving the screen horizontally and vertically.

When the black line
is on the right

When the black line
is on the left

When the black line
is at the bottom


When the black line
is at the top

Note

When the printout size is set to WIDE, the screen size cannot be adjusted vertically. (see "Setting the DIP switches" page 84)


Printer Initial Setup (continued)

- 1** When the memory image is displayed on the screen, press the SOURCE/MEMORY button.
The image from the video source appears.


- 2** Press the MENU button.
The following screen appears.

Video monitor screen


Main Menu screen

Printer window display


Part of the main menu

- 3** Select SET UP by pressing the ^ or v button.

Move the cursor to SET UP by pressing the ^ or v button.

Press the ^ or v buttons until SET UP appears.

- 4** Press the > button.
The following screen appears.

```

SET UP
▶ H SHIFT : 0 DOT
V SHIFT : 0 LINE
RM2 MODE : M&P/CYC/M&S
RM2 ACT : 1ACT/2ACT
MEM CLR : PUSH [>]
SIG ADJ : PUSH [>]
MONI ADJ : PUSH [>]

```

SET UP sub menu


```

H SHIFT :
0 DOT

```

Part of the SET UP sub menu

- 5** Select H SHIFT by pressing the ^ or v buttons, when the black line appears on the right or left.
Select V SHIFT by pressing ^ or v buttons, when the black line is at the top or bottom.

- 6** Adjust the horizontal value or vertical value by pressing the < or > button.

Item selected in step 5	The position where the black line appears	Button to be used	Operation
H SHIFT (horizontal direction)	On the right	> button	Shifting the image to the right by up to 14 dots in step 2 dots
	On the left	< button	Shifting the image to the left by up to 14 dots in step 2 dots
V SHIFT (vertical direction)	At the top	> button	In frame mode, shifting the image up by up to 6 lines in step 2 lines In field mode, shifting the image up by up to 3 lines in step 1 line
	At the bottom	< button	In frame mode, shifting the image down by up to 6 lines in step 2 lines In field mode, shifting the image down by up to 3 lines in step 1 line

- 7** Press the MENU button.
The regular screen appears.

Printer Initial Setup (continued)

To check the adjustment result

Any black line is also stored in memory with the previous image. Thus, store a new image to the memory and print it to check whether the black line disappears.

Note


When a black line still remains even after adjusting H SHIFT or V SHIFT, change the screen size. (see "Setting the printout size" page 85)

Compensating for the Input Signals

Video image recorded under poor conditions may be of poor color quality. If the signal is an PAL composite video signal or separate luminance (Y) and chrominance (C) signals, you can correct the color of the input signal to a certain extent.

Note


You cannot adjust the image once it has been stored in memory. Make adjustments before storing the image or store an image after adjustment.

- 1** When the memory image is displayed on the screen, press the SOURCE/MEMORY button.
The image from the video source appears.


- 2** Press the MENU button.
The following screen appears.

Video monitor screen


Main Menu screen

Printer window display


Part of the main menu

- 3** Select SET UP by pressing the ^ or v button.


Move the cursor to SET UP by pressing the ^ or v button.


Press the ^ or v buttons until SET UP appears.

- 4** Press the > button.
The following screen appears.


SET UP sub menu


Part of the SET UP sub menu

- 5** Select SIG ADJ by pressing the ^ or v button.

Move the cursor to SIG ADJ by pressing the ^ or v button.


Press the ^ or v button until SIG ADJ appears.


Continue to next page →

Printer Initial Setup (continued)

- 6** Press the > button.
The following screen appears.

SIGNAL ADJUST sub menu

Part of the SIGNAL ADJUST sub menu

- 7** Perform the adjustments as follows.
- ① Select the adjustment item by pressing the ^ or v button.
 - ② Perform the adjustment by pressing the < or > button.

To adjust the color intensity, hue and contrast

Adjustment	Menu adjustment item	Button and operation result	
		> button	< button
Color intensity	COLOR ^{a)}	The color intensity strengthens.	The color intensity weakens.
Contrast	GAIN	The contrast strengthens.	The contrast weakens.
Brightness	OFFSET	Becomes brighter	Becomes darker.

a) Adjust the color such that skin tones appear natural.

Note

When the black and white signal is input, set COLOR to -8. If COLOR is set to the value other than -8, display on the video monitor screen may not be colored.

When the printout or monitor image appears unnaturally bright

Adjust the input signal to the optimum level for printing.

Menu adjustment item	Selection	When selecting
AGC (Automatic gain control)	ON	Normal
	OFF	When the printout or monitor image appears unnaturally bright.

- 8** Press the MENU button.
The regular screen appears.

Matching the Video Monitor Color to the Printer Color

To match the color of the monitor image to that of the printout, adjust the monitor and printer colors such that the monitor color is the same as that of the printout. The printer outputs either of two kinds of video signals according to the MONITOR setting of the SET UP menu.

- **E TO E:** Signals are output to the monitor after being processed by the printer's circuitry.
- **THRU (through):** Signals are output to the monitor as is.


At the factory, the printer is adjusted such that the images for both signals appear identical. If they appear different, the video monitor may not be set correctly, even though the printer controls are set correctly. They may have been adjusted for another input signal. Adjust them to suit the new input signal.

First, adjust the monitor's color controls, then store and print an image. Compare the color on the monitor image with that of the printout.

- ① Adjust the monitor color by selecting the THRU (THROUGH) signal.
- ② Adjust the printer color by selecting the E TO E signal.


Printer Initial Setup (continued)

Adjusting the monitor color


- 1 Press the MENU button.
The following screen appears.

Video monitor screen


Main Menu screen

Printer window display


Part of the main menu

- 2 Select SET UP by pressing the ^ or v button.


Move the cursor to SET UP by pressing the ^ or v button.


Press the ^ or v buttons until SET UP appears.

- 3 Press the > button.
The following screen appears.

SET UP sub menu

Part of the SET UP sub menu

4 Select MONI ADJ by pressing the \wedge or \vee button.

```
MEM CLR : PUSH [>]
SIG ADJ : PUSH [>]
MONI ADJ : PUSH [>]
```

Move the cursor to MONI ADJ by pressing the \wedge or \vee button.

```
MONI ADJ :
PUSH [>]
```

Press the \wedge or \vee button until MONI ADJ appears.

5 Press the $>$ button.
The following screen appears.

```
MONITOR ADJUST
▶MONI RED : 0.....
MONI BLUE : 0.....
CHECK : PUSH [>]
DISPLAY : OFF/ON
MONITOR : E TO E/THRU
SYNC ON G : OFF/ON
```

MONITOR ADJUST sub menu

```
MONI RED : 0
----- *
```

Part of the MONITOR ADJUST sub menu

6 Select MONITOR by pressing the \wedge or \vee button.

```
MONITOR ADJUST
MONI RED : 0.....
MONI BLUE : 0.....
DISPLAY : OFF/ON
▶MONITOR : E TO E/THRU
SYNC ON G : OFF/ON
```

Move the cursor to MONITOR by pressing the \wedge or \vee button.

```
MONITOR :
<E TO E>THRU
```

Press the \wedge or \vee button until MONITOR appears.

7 Select THRU by pressing the $<$ or $>$ button.

```
MONITOR ADJUST
MONI RED : 0.....
MONI BLUE : 0.....
DISPLAY : OFF/ON
▶MONITOR : E TO E/THRU
SYNC ON G : OFF/ON
```

Switch to green.

```
MONITOR :
<E TO E>THRU
```

8 Press the MENU button.
The regular screen appears.

The image of the signal directly from the signal source, which does not pass through the printer circuit, is displayed on the video monitor.

Note


When screen display appears on the video monitor, the memory image is displayed on the monitor. Display the image from the video equipment on the video monitor by pressing the SOURCE/MEMORY button.

9 Adjust the color of the video monitor by using the monitor controls.

Printer Initial Setup (continued)


Adjusting the color of the printer's output signal

After adjusting the color of the video monitor, adjust the color of the printer's output signal.

- 1** Set MONITOR ADJ to E TO E by performing steps 1 to 7 in “Adjusting the monitor color” on page 70.
- 2** Store a new image into memory and make a printout.
You will adjust the settings by comparing the printout with the image on the monitor.
- 3** Press the MENU button.
The following screen appears.

Video monitor screen

Main Menu screen


Printer window display

PRINT QTY :

1

Part of the main menu

- 4** Select SET UP by pressing the \wedge or \vee button.

Move the cursor to SET UP by pressing the \wedge or \vee button.

SET UP :
PUSH [\wedge]

Press the \wedge or \vee buttons until SET UP appears.

- 5** Press the > button.
The following screen appears.

```
SET UP
▶H SHIFT : 0 DOT
V SHIFT : 0 LINE
RM2 MODE : M&P/CYC/M&S
RM2 ACT : 1ACT/2ACT
MEM CLR : PUSH [>]
SIG ADJ : PUSH [>]
MONI ADJ : PUSH [>]
```

SET UP sub menu

```
H SHIFT : 0 DOT
```

Part of the SET UP sub menu

- 6** Select MONI ADJ by pressing the ^ or v button.

```
MEM CLR : PUSH [v]
SIG ADJ : PUSH [v]
▶MONI ADJ : PUSH [v]
```

Move the cursor to MONI ADJ by pressing the ^ or v button.

```
MONI ADJ : PUSH [>]
```

Press the ^ or v button until MONI ADJ appears.

- 7** Press the > button.
The following screen appears.

```
MONITOR ADJUST
▶MONI RED : 0.....
MONI BLUE : 0.....
CHECK : PUSH [>]
DISPLAY : OFF/ON
MONITOR : E TO E/THRU.
SYNC ON G : OFF/ON
```

MONITOR ADJUST sub menu

```
MONI RED : 0
.....*
```


Part of the MONITOR ADJUST sub menu

Continue to next page →

Printer Initial Setup (continued)

8 Adjust the monitor while comparing the printout with the image on the monitor.

- ① Select the adjustment item by pressing the \wedge or \vee button.
- ② Perform the adjustment by pressing the $<$ or $>$ button.

Adjustment	Menu adjustment item	Button and operation result	
		$>$ button	$<$ button
Adjustment of the red component	MONI RED	The red intensity increases.	The red intensity decreases.
Adjustment of the blue component	MONI BLUE	The blue intensity increases.	The blue intensity decreases.


9 Press the MENU button.
The regular screen appears.

Note

This adjustment is performed to adjust the color of the printer output signal when the monitor color is unsatisfactory. This adjustment does not affect the printout itself. To adjust the color of the printouts, see "Adjusting the Printout Quality" on page 57.


Erasing the Screen Display

You can erase a screen display with the menu, when, for example, it is hard to see the image that is hidden behind the screen display (C, QTY, VIDEO, and others). The printer operation is identical, regardless of whether screen display are displayed on the screen. You can always see same messages as the screen display on the printer display.


- 1** Press the MENU button.
The following screen appears.

Video monitor screen


Main Menu screen

Printer window display


Part of the main menu

- 2** Select SET UP by pressing the \wedge or \vee button.

Move the cursor to SET UP by pressing the \wedge or \vee button.

Press the \wedge or \vee buttons until SET UP appears.

Printer Initial Setup (continued)

- 3** Press the > button.
The following screen appears.

```

SET UP
▶ H SHIFT : 0 DOT
V SHIFT : 0 LINE
RM2 MODE : M&P/CYG/M&S
RM2 ACT : 1ACT/2ACT
MEM CLR : PUSH [>]
SIG ADJ : PUSH [>]
MONI ADJ : PUSH [>]
  
```

SET UP sub menu

```

H SHIFT :
                                0 DOT
  
```

Part of the SET UP sub menu

- 4** Select MONI ADJ by pressing the ^ or v button.

```

MEM CLR : PUSH [>]
SIG ADJ : PUSH [>]
MONI ADJ : PUSH [>]
  
```

Move the cursor to MONI ADJ by pressing the ^ or v button.

```

MONI ADJ :
                                PUSH [>]
  
```

Press the ^ or v button until MONI ADJ appears.

- 5** Press the > button.
The following screen appears.

```

MONITOR ADJUST
▶ MONI RED : 0.....
MONI BLUE : 0.....
CHECK : PUSH [>]
DISPLAY : OFF/ON
MONITOR : E TO E/THRU
SYNC ON G : OFF/ON
  
```

MONITOR ADJUST sub menu

```

MONI RED : 0
----- * -----
  
```

Part of the MONITOR ADJUST sub menu

- 6** Select DISPLAY by pressing the ^ or v button.

```

MONITOR ADJUST
MONI RED : 0.....
MONI BLUE : 0.....
CHECK : PUSH [>]
DISPLAY : OFF/ON
MONITOR : E TO E/THRU
  
```

Move the cursor to DISPLAY by pressing the ^ or v button.

```

DISPLAY :
                                OFF<ON>
  
```

Press the ^ or v button until DISPLAY appears.

7 Select OFF by pressing the < or > button.

```
MONITOR ADJUST
MONI RED : 0.....
MONI BLUE : 0.....
CHECK : PUSH [>]
▶ DISPLAY : OFF/ON
MONITOR : E/TO E/THRU
```

Switch to green.

DISPLAY :
<OFF>ON

To display screen message

In step 7, select ON.

Note

If you set the printer output signal specification to THRU (through), screen display does not appear, even when you switch ON to green.

8 Press the MENU button.
The regular screen appears.

Printer Initial Setup (continued)


Selecting the Operation Mode for Automatic Printing Capabilities

You can control the printer with the RM-91 remote commander connected to the REMOTE 2 connector on the rear panel.

In addition to the above, the printer can be controlled by the pulse signal input to REMOTE 2. (see page 92)


Note

When controlling the printer with the remote commander (not supplied), set INTERVAL of the memory mode to OFF. (see page 27)

- 1 Press the MENU button.
The following screen appears.

Video monitor screen

Main Menu screen


Printer window display

PRINT QTY :

1

Part of the main menu

- 2** Select SET UP by pressing the \wedge or \vee button.


Move the cursor to SET UP by pressing the \wedge or \vee button.

SET UP :
PUSH [\vee]

Press the \wedge or \vee buttons until SET UP appears.

- 3** Press the $\>$ button.
The following screen appears.


SET UP sub menu

H SHIFT :
0 DOT

Part of the SET UP sub menu

- 4** Select RM2 MODE by pressing the \wedge or \vee button.

Move the cursor to RM2 MODE by pressing the \wedge or \vee button.



RM2 MODE :
<M&P> CYC/M&S

Press the \wedge or \vee button until RM2 MODE appears.

Continue to next page →

Printer Initial Setup (continued)

5 Select the desired operation method by pressing the < or > button.



Position < > over the desired operation method.

Switch the desired operation method to green.

Type of control operation	Operation method
M & P (MEMORY & PRINT)	Storing an image into memory page and printing memory image. If you press the switch of the remote commander while printing is being performed, the selected image is queued, with the queued image being printed once printing has been completed.
CYC (CYCLIC MEMORY)	Storing images to memory page cyclically whenever you press the switch of the remote commander. The printer continues to store images, replacing previously stored images with the new one.
M & S (MEMORY & STOP)	Storing an image to memory page whenever you press the switch of the remote commander. The printer stops storing images to memory page once images have been stored to all memory pages. The Message STOP STOP STOP appears.



6 Select RM2 ACT by pressing the ^ or v button.

Press the ^ or v button until RM2 ACT appears.

Move the cursor to RM2 ACT by pressing the ^ or v button.

7 Select the desired operating condition by pressing the < or > button.

Position < > over the desired operating condition.

Switch the desired operating condition to green.

Type of operating condition	Operating condition
1ACT	Whenever you press the switch of the remote commander, the printer automatically stores the image into the next memory page. You cannot check an image to be next stored in memory.
2ACT	Whenever you press the switch of the remote commander, the printer automatically stores the image and you can check an image to be next stored in memory.

8 Press the MENU button.
The regular screen appears.

Printer Initial Setup (continued)

Controlling the printer by the remote commander

The printer is remotely controlled according to the RM2 ACT setting.

When ACT1 is selected

Operation method selected in step 5	Operation
M & P (MEMORY & PRINT)	Storing an image in free memory page and simultaneously printing it whenever you press the switch of the remote commander. If you press the switch of the remote commander while printing is being performed, the selected stored image is queued, with the queued image being printed, once printing has been completed.
CYC (CYCLIC MEMORY)	Storing an image to free memory page cyclically whenever you press the switch of the remote commander. The printer continues to store images, replacing previously stored images with the new one.
M & S (MEMORY & STOP)	Storing an image to free memory page whenever you press the switch of the remote commander. The printer stops storing images to memory page once images have been stored to all memory pages. The message STOP STOP STOP appears.

To make the message STOP STOP STOP disappear

Press the STOP button. The printer is reset to the normal printing mode.

When SPLIT MEM of the MULTI PICTURE sub menu is set to 4 or 16

For example, when the memory mode is set to store four reduced images in a memory page, press the switch of the remote commander to store one reduced image to each position. When M & P is selected, the printer stores the image for the fourth position and make a printout of four reduced images when you press the switch of the remote commander the fourth time.

Note

When the image in a memory page is queued for printing, that memory page is automatically skipped.

When ACT2 is selected

Display the video source image before starting operation.

Operation method selected in step 5

M & P (MEMORY & PRINT)

When SPLIT MEM is set to OFF

- ① Press the switch of the remote commander at the instant when the image you want to print appears.
The image is stored in the free memory page.
- ② Press the switch of the remote commander.
The image stored in step 1 is printed.
The video source image appears on the video monitor.
- ③ Repeat steps 1 and 2.

When SPLIT MEM is set to 4 or 16

Once the printer have been stored the last reduced image, the printer makes a printout of multiple reduced images. For example, when 4 is selected, the fourth reduced image is stored after you press the switch of the remote commander the seventh time. When you press the switch of the remote commander the eighth time, the printer makes a printout of the four reduced images.

CYC (CYCLIC MEMORY)^{a)}

- ① Press the switch of the remote commander at the instant when the image you want to print appears.
The image is stored in the free memory page.
- ② Press the switch of the remote commander.
The video source image appears on the video monitor.
- ③ Repeat steps 1 and 2.
An image is stored to the memory page cyclically. The printer continues to store images, replacing previously stored images with the new one.

M & S (MEMORY & STOP)^{a)}

- ① Press the switch of the remote commander at the instant when the image you want to print appears.
The image is stored in the free memory page.
 - ② Press the switch of the remote commander.
The video source image appears on the video monitor.
 - ③ Repeat steps 1 and 2.
An image is stored to the memory page whenever you press the switch of the remote commander. The printer stops storing images to the memory page once images have been stored to all memory pages. The message STOP STOP STOP appears.
-

a) When SPLIT MEM is set to 4 or 16: The printer repeats the operation described above, whenever the printer stores one reduced image to each position.

To make the message STOP STOP STOP disappear

Press the STOP button. The printer is reset to the normal printing mode.

Note


When the image in a memory page is queued for printing, that memory page is automatically skipped.

Printer Initial Setup (continued)

Setting the DIP Switches

You can make the following five settings by using the DIP switches on the rear panel.

- Setting of the printout size (see page 85)
- Setting of the baud rate for computer communication (see page 85)
- Setting of the MENU button function ON/OFF (see page 85)
- Selecting whether lock to an external synchronous signal (GENLOCK) (see page 86)
- Selecting whether the operation and error tones sound (see page 86)

Note

Before changing the DIP switch setting, turn the printer power off. The former setting with the DIP switch remains if you change the switch settings with the power on.

Setting the printout size

When you print an image that is narrower or wider than the standard screen size, the black frame may be printed or the image may be partially cut. In such a case, you can change the screen size.

The printer supports the following three sizes, NA (NARROW), NO (NORMAL) and W (WIDE).

Use DIP switches SW1 and SW2 to set the printout size.

When changing	Printout size	Size (dots × line)	SW1 set position	SW2 set position
When a black frame is printed.	NA (NARROW)	708 (H) × 448 (V)	OFF	ON
Normal	NO (NORMAL)	720 (H) × 472 (V)	OFF	OFF
When an image is partially cut	W (WIDE)	772 (H) × 488 (V)	ON	OFF

Setting the baud rate for computer communication

When controlling the printer with a personal computer connected to the RS-232C connector, select the appropriate baud rate. Use DIP switches SW3 and SW4 to set the baud rate.

Baud rate (bps)	SW3 set position	SW4 set position
1200	OFF	OFF
2400	ON	OFF
4800	OFF	ON
9600	ON	ON

Selecting whether to enable the MENU button functions

You can inhibit the MENU button operation. To make printouts with fixed menu values without performing any menu operations during normal printing, set the MENU button function to OFF. Use DIP switch SW5 to select the MENU button function.

When changing	MENU button function	SW5 set position
To make the MENU button operable	OFF (NORMAL)	OFF
To inhibit MENU button operation	ON (INHIBIT)	ON

Printer Initial Setup (continued)

Selecting whether to lock to an external synchronous signal (GENLOCK)

You can select whether the printer is locked to the external synchronous signal when sending a video signal to the video monitor connected to the OUTPUT connectors on the rear panel.

Use DIP switch SW6 to set this item.

When changing	GENLOCK ON/OFF	SW6 set position
Not to be locked to the external synchronous signal	OFF (FREE)	OFF
To be locked to the external synchronous signal	ON (LOCK)	ON

Note

When GENLOCK is not to ON (LOCK), the printer may not work correctly if irregular signal (such as search mode or stop mode signal) is input to the printer.

Selecting whether the operation and error tones sound

Whenever you use the MEMORY IN, SOURCE/MEMORY, PRINT and STOP button on the front panel, the operation tone sounds (one time). When an error occurs, the alarm tone sounds (three times). You can turn these tones on and off. Use DIP switch SW7 to set this item.

When changing	BEEP (tone) ON/OFF	SW7 set position
To sound the tones	ON (sounding)	OFF
To disable the tones	OFF (not sounding)	ON


Note

Even when BEEP is set to OFF (not sounding), operation tone or error tones sound when the printer power is turned on, or any error has occurred.


Displaying the DIP switch settings on the video monitor

You can check the DIP switch settings on the video monitor and printer window display.

The following operation can only be performed by using the buttons on the printer.

Press the SOURCE/MEMORY button holding down the STOP button.
The following screen appears.

The selected setting is displayed in green.

To return to the regular screen, press the MENU button.

Note

When displaying the DIP switch settings on the video monitor and in the printer window display:

Characters and numbers displayed on the video monitor display and in the printer window display are abbreviated. Note that some of the DIP switch settings on the rear panel of the printer differ from those displayed on the video monitor display and in the printer window display or those used for explanations of the DIP switch settings on pages 84 to 86 of this manual.

See the following table.

Indication on monitor display and in the printer window display		DIP switch setting on the rear panel
MENU	ALLOW	OFF (NORMAL)
	MASK	ON (INHIBIT)
BEEP	BEEP	ON (sounding)
	MASK	OFF (not sounding)

Precautions

Safety

- This unit requires 220 to 240 V AC power supply.
- Be careful not to damage the power cable by placing or dropping heavy objects on it; it is dangerous to use the unit with a damaged cable.
- If you do not intend to use the unit for a long time, disconnect the power cable.
- Unplug the power cable by holding the plug, not the cable itself.
- Do not disassemble the unit.
- Do not remove the cover. There is a danger of electric shock from the internal parts.
- Be careful not to spill water or other liquids on the unit, or to allow combustible or metallic material to enter the cabinet. If used with foreign matter in the cabinet, the unit is liable to fail, or be a cause of fire or electric shock.
- Vents are provided to prevent the unit from overheating. Be careful not to obstruct them with other units or by covering the unit with a cloth etc.
- If the unit should malfunction or if a foreign body falls into the cabinet, disconnect the power immediately and consult your Sony dealer.

Installation

Avoid placing the unit in a location subject to:

- mechanical vibration
- high humidity
- excessive dust
- direct sunlight
- extremely high or low temperature


Do not use either electronic equipment near the unit. The unit will not work properly in strong electromagnetic fields.

On condensation

- If the printer is subjected to wide and sudden changes in temperature, such as when it is moved from a cold room to a warm room or when it is left in a steamy or damp room due to the heater, condensation may form inside the printer. This will cause its parts to damage and spots to appear on the printout. Therefore, avoid using the printer in locations subject to wide and sudden changes of temperature.
- If the printing pack (ink ribbon cassette and paper) is subjected to wide and sudden changes in temperature, such as when it is moved from a cold room to a warm room, condensation may form on the ink ribbon or print paper. This will cause the printer to malfunction. Also, if the printing pack is used in this state, spots are likely to appear on printout. Therefore avoid the printing pack in locations subject to wide and sudden changes of temperature.
- To store a half-used printing pack, replace it in this original packing and reseal the package. If possible, keep the sealed printing pack in a cool, dark location. To subsequently use the printing pack, place it in the sealed package, in a warm room for several hours. Doing so prevents condensation from forming when the printing pack is removed from its package.

Location

To prevent internal heat built-up, leave enough room around the printer for air to circulate through the vents on the left.

Cleaning

Clean the cabinet, panel and controls with a dry soft cloth, or a soft cloth lightly moistened with a mild detergent solution. Do not use any type of solvent, such as alcohol or benzene, which may damage the finish.

Ink Ribbon Cassette and Paper

The video printer is designed to use ink ribbon cassette and paper (not supplied) with VP mark to obtain a high printing quality. The ink ribbon cassette and paper are sold together as sets; a set for color printing on paper, for monochrome printing on paper, and for monochrome printing on laminated paper. Use whichever cassette/paper set is appropriate:

Color printing pack UPC-1010 (not supplied)

One pack contains an ink ribbon cassette and 100 sheets of $4\frac{1}{4} \times 5\frac{3}{4}$ " (A6 size) paper.

Monochrome printing pack UPC-1020 (not supplied)

One pack contains a B & W ink ribbon cassette and 100 sheets of $4\frac{1}{4} \times 5\frac{3}{4}$ " (A6 size) paper.

Self-laminating color printing pack UPC-1040 (not supplied)

One pack contains an ink ribbon cassette and 75 sheets of $4\frac{1}{4} \times 5\frac{3}{4}$ " (A6 size) paper for automatic laminate coating.

Specifications

Power requirements

220 to 240 V AC, 50/60 Hz

Power consumption

1.0 A max.

at 25°C, 240 V AC (~)

Operating temperature

5°C to 40°C (41°F to 104°F)

Operating humidity

20% to 80% (no condensation allowed)

Storage and transport temperature

-20°C to 60°C (-4°F to 140°F)

Storage and transport humidity

20% to 90% (no condensation allowed)

Dimensions

About 424 × 91 × 397 mm

(w/h/d) (16³/₄ × 3⁵/₈ × 15³/₄ inches)

Mass

About 8.5 kg (18 lb 12 oz)

Printing system

Sublimation heat transfer printing

Thermal head

6.72 dot/mm (608 dots)

Total gradation

256 levels each for yellow, magenta, and cyan

Frame memory

UP-1800EPM: One frame memory

UP-1850EPM: Three frame memories

Printing time

Approximately 60 seconds (normal size color printing)

Approximately 30 seconds (monochrome printing)

TV system

PAL B.G.I. standards

Input connectors

RGB SYNC

(analog RGB signal): BNC connector
0.7 Vp-p, 75 ohms (75 ohm termination switch set to ON)

S-VIDEO (Separate luminance (Y) and chrominance (C) signals): 4-pin mini-DIN

Y: 1 Vp-p

C: 0.3 Vp-p color burst

75 ohms (75 ohm termination switch set to ON)

VIDEO (PAL composite video signal): BNC connector

1 Vp-p, 75 ohms (75 ohm termination switch set to ON), sync negative

AC IN (for power input)

Output connectors

RGB SYNC (analog RGB signal): BNC connector

RGB: 0.7 Vp-p, 75 ohms
(75 ohm termination switch set to ON)

SYNC: 1 Vp-p

S-VIDEO (Separate luminance (Y) and chrominance (C) signals):

4-pin mini-DIN

Y: 1 Vp-p, 75 ohms

C: 0.3 Vp-p Color burst, 75 ohms (75 ohm termination switch set to ON)

VIDEO (PAL composite video signal): BNC connector

1 Vp-p, 75 ohms (75 ohm termination switch set to ON), sync negative

Controls connectors

REMOTE 1 (front panel, for the supplied remote commander only): special mini jack

REMOTE 2 (automatic printing connector, for the remote commander (not supplied)):

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(see "Using the automatic printing capabilities" page 92)

RS-232C (Computer control interface):

D-SUB 25-pin connector

Ink ribbon cassette and paper sets

Color printing pack: UPC-1010 (100 sheets)

B & W printing pack:

UPC-1020 (100 sheets)

Self laminating color printing pack:

UPC-1040 (75 sheets)

Supplied accessories

Color printing pack UPC-1010 (1)

Paper tray (1)

Paper cover (1)

Remote control unit RM-5100 (1)

Connecting cable for the remote control unit (1)

Dry battery SUM-3 (NU) (2)

AC power cord (1)

Instructions For Use (1)

Specifications (continued)


Using the automatic printing capabilities (REMOTE 2)

If you send the $\overline{\text{MIN}}$ remote control pulse signals illustrated below through the REMOTE 2 connector, the printer is remotely controlled according to the remote control setting. (see page 78).


Turn on the power of the printer and select the input signal. Set the monitor display to the input signal. Send a remote control pulse signal at timing shown below:

Regulations of remote control pulse

MEMORY IN TIMING

PRINT TIMING

Design and specifications are subject to change without notice.

Error Messages

If a problem occurs, the ALARM lamp lights in orange when an error message or a warning message appears. The message which states the problem appears on the monitor and in the printer window display. This section lists error and warning messages in alphabetical order, together with their possible causes and remedies. Note the message and act accordingly.

Error message	Possible causes and remedies
CHECK RIBBON SETTING	The front panel (on the right from the user's standpoint) opens accidentally during printing. → Close the front panel. (see page 9)
FEED ERROR	The paper jams as it is being fed into the ribbon area around the paper tray. → Remove the jammed paper from the printer. (see page 94)
HEAD IN COOLING	The thermal head has overheated. → Leave the printer idle and until the head cools and this message disappears.
NO CARTRIDGE	The ink ribbon cassette is not correctly installed. → Insert the ink ribbon cassette correctly. (see page 8)
NO PAPER	The paper has been exhausted. → Load paper. (see page 10)
PAPER ERROR	The paper that cannot be used with this printer has been loaded. → Insert the appropriate paper. (see page 90)
PREHEATING	The thermal head is preheating. → Leave the printer until the head has preheated and this message disappears.
REMOVE PRINTS	The paper has jammed near the paper cover. → Remove the jammed paper from the printer. (see page 94)
REMOVE STUCK PAPER	The paper has jammed during printing. → Remove the jammed paper from the printer. (see page 94)
RIBBON & PAPER MISMATCH	The ink ribbon cassette and paper are not compatible. → Use a compatible cassette/paper combination. (see page 90)
RIBBON DOOR OPEN	The front panel (on the right from the user's standpoint) is open. → Close the front panel. (see page 9)
RIBBON END	The ink ribbon cassette has been exhausted. → Insert a new ink ribbon cassette. (The ink ribbon cannot be reused.) (see page 8)
RIBBON ERROR	An ink ribbon cassette that cannot be used with this printer has been loaded. → Insert the appropriate ink ribbon cassette. (see page 90)

Error Messages (continued)

If ERRORxx appears

If this message appears, perform the following.

- 1** Turn off the power of the printer.
- 2** Remove the ink ribbon cassette, paper cover and paper tray, and check for any paper jams inside the printer. (see “Loading an Ink Ribbon Cassette” page 8 and “Loading Paper” page 10)
If you find any jammed paper, remove it carefully.

If the ink ribbon cassette cannot be removed, or the jammed paper cannot be removed, contact your Sony service facility.

- 3** Insert the ink ribbon cassette, paper cover and paper tray to the printer.
- 4** Turn on the power of the printer.
When the message does not appear, you can use the printer as normal.
However, the image stored to memory will have been cleared. Store the image to memory again.

If the same message appears again, the printer must not be operated. Turn off the power immediately and contact your Sony service facility.


If the paper jams

If the paper jams as it is being fed into the ribbon area during printing, or when being fed into the paper cover area, printing stops and a message appears on the monitor and in the printer window display, according to where the jam has occurred.


Message	Position where the paper has jammed
FEED ERROR	Before printing and being fed into the ribbon area
REMOVE STUCK PAPER	During printing, inside the printer
REMOVE PRINTS	Instantaneously before completing printing, near the paper cover

When FEED ERROR appears


- 1** Remove the paper cover.
When any printouts have been ejected on the paper cover, remove those printouts first before removing the paper cover.


- 2** Check whether any paper is stuck inside the printer. If you find a jammed sheet, slowly pull it into the paper tray.

- 3** Remove the paper tray.


- 4** Load the paper into the paper tray correctly.
Discard the paper removed in step 2.

Continue to next page →

Error Messages (continued)

- 5** Slide the paper tray and paper cover back into the printer.

When REMOVE STUCK PAPER appears

Perform the same operation as that performed when FEED ERROR appears. When you cannot remove the jammed paper, remove the ink ribbon cassette, too. If you find a jammed sheet inside the printer, remove it carefully.

When REMOVE PRINTS appears

Carefully remove the jammed paper from near the paper cover.

Troubleshooting


The following troubleshooting checks will help you correct the most common problems encountered with your unit. Before proceeding with these trouble checks, first check that the power cord is firmly connected. Should the problem persist, unplug the unit and contact your Sony dealer or local authorized Sony service facility.

Symptom	Possible causes and remedies
Nothing appears on the monitor.	<p>The POWER switch of the printer is not set to ON. → Set the POWER switch of the printer to ON.</p> <p>The POWER switch of the monitor is not set to ON. → Set the POWER switch of the monitor to ON.</p> <p>Connections may not be correct. → Make connections correctly. (see page 52)</p>
Any message does not appear on the regular screen.	<p>If an incorrect sync signal is input, nothing may appear on the monitor. → In this case, check the monitor first by pressing the SOURCE/MEMORY button to display the image stored in memory. If an image appears, the monitor is working correctly. Change the INPUT SEL settings on the menu screen. (see page 12) Or, set the connected video equipment to playback mode, if it is in another mode such as stop mode.</p>
Any message and image do not appear on the regular screen.	<p>If an image stored in memory appears when the SOURCE/MEMORY button is pressed, the MONITOR settings on the MONITOR ADJUST sub menu is set to THRU. → Change the MONITOR settings on the MONITOR ADJUST sub menu to E TO E. (see page 71)</p>
The printer does not print.	An error message appears on the display. (see page 93)
A black line appears on the printout.	<p>A portion corresponding to there being no signal is printed in black. → Shift the printout area. (see page 63) Store a new image and print it.</p>
The printer makes a printout with black frame.	<p>A portion corresponding to there being no signal is printed in black. → Change the printout size to make it narrow. (see page 85) Store a new image and print it.</p>
The printed image is partially cut out.	<p>Only a part of video signal has been stored. → Change the printout size to make it wide. (see page 85) Store a new image and print it.</p>
A caption is not printed clearly.	<p>Printed in field mode. → Store the image in frame mode and print it in frame mode.</p>
The printout is blurred.	<p>The quickly moving image has been stored in frame mode. → Change the mode to field mode, then print it again.</p>

Location and Function of Parts and Controls

For details, refer to the pages indicated in (00).

Front

1 ① POWER switch

Press to turn the printer on or off.

2 Paper tray/paper cover (54)

Paper tray: Load paper into this tray.

Paper cover: The printout is ejected to this tray

3 Printer window display (102)

This window displays the messages that also appear along the bottom edge of the monitor screen. It also displays the menu screen line to which the cursor is positioned. If an error occurs, corresponding error message is displayed.

4 MENU button

This button is used to display menus or to return to the regular screen from the main menu or sub menus.

The DIP switch setting allows you to make this MENU button disable to work.

5 EXEC button (31, 43, 47, 50)

Press this button to return to the previous menu. Also, this button is used to erase the memory image or enter characters for a caption.

6 SOURCE/MEMORY button (15, 32, 34, 45, 48, 64, 66, 71, 87)

Press to select which signal is to be output to the monitor.

The memory image and source image are changed whenever you press this button.

7 ⇨ MEMORY IN button (15, 32, 34)

Press to store an image into memory.

8 ALARM lamp (93)

This lamp lights, in orange, when the paper has jammed or another error occurs.

9 ☐ PRINT button (16, 17, 33, 34)

Press to make printouts.

10 PUSH OPEN button (8)

Press to open the right front panel door when loading an ink ribbon cassette.

11 Remote sensor (56)

Aim the head of the remote control unit toward this sensor.

12 REMOTE 1 connector (54)

Used to Connect the remote control unit (supplied) when being used as a wired type.

13 STOP button (16, 33, 82, 83, 87)

Press to stop printing midway.

Press this button when the message "STOP STOP STOP" appears.

14 MEMORY PAGE button (25, 38, 47, 50)

Press to select the memory page.

15 Cursor keys

Press to position the cursor. Select a desired item from the menu by pressing the ^ or v button and set the value by pressing the < or > button.

Also, these keys are used to enter characters for a caption.

16 LCD control


Turn to adjust the contrast of the printer window display. To increase the contrast, turn the control clockwise. To reduce the contrast weaker, turn the control counterclockwise.

17 PUSH indication (10)

Press this position to remove the paper tray.

Location and Function of Parts and Controls (continued)

Rear

1 INPUT connectors (52)

Used to connect to the video equipment for source image.

Connector	Connectable equipment
S-VIDEO	Video equipment with a Y/C separated output
VIDEO	Video equipment with a composite video signal output
RGB SYNC	Video equipment with RGB/SYNC connectors

See "Important safeguards/notices for use in the medical environment" on page 2.

2 75-ohm termination switch (for RGB input signal and composite video signal) (52)

Normally, set this switch to ON. Set it to OFF if the input signal should drop when you connect additional equipment to the video equipment.

3 REMOTE 2 connector (54)

Used to connect the RM-91 remote commander (not supplied) or input remote control pulse signals for automatic printing.

4 RS-232C connector (54)

Used to connect the computer used to control the printer.

See "Important safeguards/notices for use in the medical environment" on page 2.

5 DIP switches (84)

Set the switches according to the printout size, baud rate setting for the RS-232C connector, MENU button enable/disable setting, whether locking/non-locking to an external synchronous signal, and whether the tones are turned on or off.

6 ~ AC IN connector (52, 53, 54)

Used to connect to a wall outlet with the supplied power cord.

7 ⚡ Equipotential terminal

Used to connect to the equipotential plug to bring the various parts of a system to the same potential.


8 OUTPUT connectors (53)

Used to connect to the video monitor.

Connector	Connectable video monitor
S-VIDEO	Video monitor with a Y/C separated input
VIDEO	Video monitor with a composite video signal input
RGB SYNC	Video monitor with RGB/SYNC connectors

See "Important safeguards/notices for use in the medical environment" on page 2.

Remote Commander RM-5100

1 PRINT QTY + and - buttons (19)


Used to set the number of copies of one printout (on the regular screen).

Button	Operation
+	Increases the number of copies.
-	Reduces the number of copies.

2 MULTI PICTURE button (26)

Press to access the MULTI PICTURE sub menu directly from the any other screen.

3 Menu control keys

a MENU button

This button is used to display menus or to return to the regular screen from the main menu or sub menus.

The MENU button can be enabled/disabled by setting a DIP switch.

b Cursor keys

Press to position the cursor. Select a desired item from the menu by pressing the ^ or v button and set the value by pressing the < or > button.

Also, these keys are used to enter characters for a caption.

c EXEC button (31, 43, 47, 50)

Press this button to return to the previous menu. Also, this button is used to enter characters for a caption.

4 PRINT button (16, 17, 33, 34)

Press to make printouts.

5 MEMORY IN button (15, 32, 34)

Press to store an image into memory.

6 SOURCE/MEMORY button (15, 32, 34, 45, 48, 64, 66, 71, 87)

Press to select which signal is to be output to the monitor.

The memory image and source image are changed whenever you press this button.

7 STOP button (16, 33, 82, 83, 87)

Press to stop printing midway.

Press this button when the message "STOP STOP STOP" appears.

8 COLOR ADJUST button (59)

Press to access the COLOR ADJUST sub menu directly from any other screen.

9 MEMORY PAGE button (25, 38, 47, 50)


Press to select the memory page.

Location and Function of Parts and Controls (continued)


Comparing Printer Window Display with the Video Monitor Display

The printer window display differs from the monitor display in that the display range is narrower and the number of characters is limited. The displayed contents, however, are the same.
A brief explanation of the printer window display is given below, in comparison with the monitor screen.

Regular screen


Displays a line of the numbers, one of which blinks in white in the multiple reduced-image area. < > is positioned on the number that blinks on the monitor display.


Displays the bottom line of the monitor screen.

Main menu screen and sub menu screen


Displays the item to which the cursor is positioned on the monitor display


Displays the selections for the item to which the cursor is positioned on the monitor display < > is positioned on the selected selection.

CAPTION sub menu screen

Displays the line to which the highlighted cursor is positioned in the character display area of the monitor display (blinks in black)

Displays the line of the character entry area to which the highlighted cursor (selected character) is positioned.


Displays a selected character (blinks in black)

Monitor Display

There are two types of screen display; the regular screen display and the menu screen.

Regular screen message

When you first turn on the printer, the regular screen message appears.

1 C (Caption)

C is displayed in white when the printer is set to print a caption consisting of the date and/or comments.

C is displayed in dark blue when the printer is not set to print a caption consisting of the date and/or comments.

M is displayed in white when the printer is set to print a mirror caption.

2 Error message display area

Error messages are displayed.

3 Warning message display area

Warning messages are displayed.

4 Number of four or 16 reduced image area

When the printer is set to store multiple reduced images into memory, corresponding numbers appear to indicate the memory status.

5 I1

Displayed when the printer is set to automatically store images at preset intervals. The number corresponds to the selection made with the INTERVAL item.

6 Image type display

This indicates the type of image shown on the monitor screen.

When the image being played back from print source equipment is displayed on the screen, the corresponding print source (the input signal connector name, for example VIDEO) appears. When an image stored in memory is displayed on the screen, MEMORY appears.

Location and Function of Parts and Controls (continued)

7 Print mode display

This indicates the selected print mode.
Several examples are shown below:

Display	Print mode
NORM	Makes a printout of one normal image
N2	Makes a printout of two identical normal images
N1_M	Makes a printout of one normal image whose printout area size is reduced to 70%.
MIR	Makes a printout of one mirror image
M16	Makes a printout of 16 reduced mirror images

8 Memory page display

The memory page you select appears.
The memory page whose image is being printed blinks. A memory page whose image is queued for printing lights in green.
The following shows several examples.

Display	Meaning
1/3FRM (UP-1850EPM) 1/1FRM (UP-1800EPM)	The first page is selected in frame mode (UP-1850EPM). The frame mode is selected (UP-1800EPM).
2/6FLD (UP-1850EPM) 2/2FLD (UP-1800EPM)	The second page is selected in field mode.

9 Number of copies to be printed


Indicates the number of copies to be printed.
This item blinks while the printer is busy. Also, the color changes to indicate the progress while making a color printout, as follows:
Printing start - yellow - magenta - cyan - printing end.
When making black and white printouts, this blinks in white.

Menu screen

Menu screen tree-chart

The menu screen configuration is shown using the tree-chart.

For details, refer to the pages indicated in ().


Note

In addition to the above menu screens, there is the DIP switch setting display. This display, however, is not used for menu operation. (see page 87)

Location and Function of Parts and Controls (continued)

Menu screen display

About remaining amount of ribbons

- Use the remaining amount of ribbons as a guide. The printer may not display remaining amount of ribbons depending on the type of ribbons.
- The printer starts counting from 25 when remaining amount of ribbons is less than 36.

Display color

The color indicates the printer status.

Display color	Meaning
Light blue	Indicates the menu name.
Green	In the menu item column, indicates the selected item. In the selection column, indicates an item that has already been set or one that must be set.
White	In both the item and selection column, indicates that the item has not been selected or has not yet been set.
Dark blue	Indicates that this item or selection cannot be selected. They are functions which become effective depending on another item or selection settings.

Index

A

- About this manual 5
- Adjustment
 - input signal 66
 - matching the video monitor color to the printer color 69
 - printout area 63
 - printout color 58
 - printout size 85
 - sharpness of the printout 57

C

- Caption
 - CAPTION sub menu 40
 - entering 41
 - making a printout with a caption 44
- Connections
 - for enabling remote control 54
 - for storing video images 52
 - for viewing images to be printed on the video monitor 53

D

- Deleting the memory image
 - all images from a single memory page 45
 - part of multiple reduced images from a signal memory page 48
- Display
 - erasing screen display 75
 - monitor display 103
 - printer window display 102

E

- Error messages 93

F

- Field mode and frame mode 22
- Full-size printout
 - making 14
 - setting 27

I

- Initial setup
 - DIP switch settings
 - baud rate for computer communication 85
 - MENU button function 85
 - printout size 85
 - selecting whether to enable the operation and error tone 86
 - selecting whether to lock to an external synchronous signal 86

setup with menu

- changing printout area 63
- compensating for input signals 66
- erasing the screen display 75
- matching the video monitor color to printer color 69
- operation mode of automatic printing 78
- Ink ribbon cassette
 - ink ribbon cassette/paper set 90
 - loading 8
 - removing 8
- Input signal
 - compensating 66
 - connection 52
 - selecting 12
- Installation
 - assembling the supplied paper cover and paper tray 54
 - connections 52
 - precautions 88
 - supplied accessory 51

J

- Jammed paper
 - error message 94
 - removing 94

L

- Location and function of parts and controls
 - front 98
 - rear 100
 - remote control unit (supplied) 101

M

- Making a printout
 - full size-printout 14
 - multiple copies of identical image 18
 - operation flowchart for making a printout 21
 - printout of multiple images from multiple memory pages 35
 - printout of multiple reduced images 31
 - printout with a caption 44
 - variations of printouts 20

Memory

- about memory 22
- deleting memory image 45
- field mode and frame mode 22
- selecting field mode and frame mode 23
- selecting the memory mode 26
- selecting a memory page 25

Memory page 24

Menu

- CAPTION sub menu 40
- COLOR ADJUST sub menu 60
- main menu 105
- MEMORY CLEAR sub menu 46, 49, 105
- MONITOR ADJUST sub menu 53, 71, 73, 76, 105
- MULTI PICTURE sub menu 27, 29, 36, 105
- PICK IMAGE sub menu 37, 105
- SET UP sub menu 65, 67, 70, 73, 76, 105
- SIGNAL ADJUST sub menu 68, 105
- tree-chart 105

Mirror-image printout 30

Multiple reduced-image printout

- making 31
- setting 27

O

On-screen menu

- how to read 106
- tree-chart 105

P

Paper

- ink ribbon cassette/paper set 90
- loading 10
- removing the jammed paper 94

Precautions

- cleaning 89
- condensation 89
- installation 88
- safety 88

Preparation

- adjustment of the printout quality 57
- confirming the supplied accessory 51
- connections 52
- DIP switch setting 84
- loading an ink ribbon cassette 8
- loading paper 10
- printer initial setup 63
- remote control unit 55

Printer

- initial setup 63
- overview 6
- system configuration 7
- printer window display 102

Printout

- composed image 35
- imperfect printout
 - with black line 63
 - with black frame 85
 - partially cut 85
 - with poor picture quality 57
- with an insert 33
- multiple reduced images 31
- with white borders 30
- without white borders 30
- mirror 30
- image size in a single printout 30

R

Remote control unit (not supplied)

- foot switch 56
- pulse signal 92
- selecting operation mode 78

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